

Nur für Forschungszwecke

EMX2-Specific Polyklonaler Antikörper

Katalog-Nr.: 19941-1-AP



Allgemeine Informationen

Katalog-Nr.:

19941-1-AP

Größe:

150ul, Konzentration: 140 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_004098

GeneID (NCBI):

2018

Vollständiger Name:

empty spiracles homeobox 2

Berechnete Masse:

28 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

EMX2 belongs to the EMX homeobox family. Transcription factor, which in cooperation with EMX2, acts to generate the boundary between the roof and archipallium in the developing brain. EMX2 may function in combinations with OTX1/2 to specify cell fates in the developing central nervous system. Defects in EMX2 are the cause of schizencephaly which is an extremely rare human congenital disorder characterized by a full-thickness cleft within the cerebral hemispheres. This antibody is specific to EMX2.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

