

HTR2C Polyklonaler Antikörper

Katalog-Nr.: 19959-1-AP

Allgemeine Informationen

Katalog-Nr.:

19959-1-AP

Größe:

150ul, Konzentration: 133 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_000868

GeneID (NCBI):

3358

Vollständiger Name:

5-hydroxytryptamine (serotonin)
receptor 2C

Berechnete Masse:

52 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus

Hintergrundinformationen

HTR2C, also named as HTR1C, belongs to the G-protein coupled receptor 1 family. HTR2C is one of the several different receptors for 5-hydroxytryptamine (serotonin), a biogenic hormone that functions as a neurotransmitter, a hormone, and a mitogen. It mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system. It is glycosylated at N-term. The antibody recognizes the C-term of HTR2C.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

