

Nur für Forschungszwecke

KCND2-Specific Polyklonaler Antikörper

Katalog-Nr.:19964-1-AP



Allgemeine Informationen

Katalog-Nr.:

19964-1-AP

Größe:

150ul, Konzentration: 160 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_012281

GeneID (NCBI):

3751

Vollständiger Name:

potassium voltage-gated channel, Shal-related subfamily, member 2

Berechnete Masse:

71 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

KCND2, also named as KIAA1044 and KV4.2, belongs to the potassium channel family and D (Shal) subfamily. It is a pore-forming (alpha) subunit of voltage-gated rapidly inactivating A-type potassium channels. KCND2 may contribute to I(To) current in heart and I(Sa) current in neurons. Channel properties are modulated by interactions with other alpha subunits and with regulatory subunits. This antibody is specific to KCND2.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

