

Nur für Forschungszwecke

transcription termination factor- Specific Polyklonaler Antikörper



Katalog-Nr.:20005-1-AP

Allgemeine Informationen

Katalog-Nr.:
20005-1-AP

Größe:
150ul

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_007344

GeneID (NCBI):
7270

Vollständiger Name:
transcription termination factor, RNA
polymerase I

Berechnete Masse:
103 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

TTF1 is a multifunctional nucleolar protein that terminates ribosomal gene transcription, mediates replication fork arrest and regulates RNA polymerase I transcription on chromatin. It plays a dual role in rDNA regulation, being involved in both activation and silencing of rDNA transcription. Interaction with BAZ2A/TIP5, TTF1 recovers DNA-binding activity. The antibody is specific to TTF1. This TTF1 is different to another TTF1 (Thyroid transcription factor 1).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

