

RASSF7 Polyklonaler Antikörper

Katalog-Nr.: 20008-1-AP

Allgemeine Informationen

Katalog-Nr.:

20008-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_003475

GeneID (NCBI):

8045

Vollständiger Name:

Ras association (RalGDS/AF-6)

domain family (N-terminal) member

7

Berechnete Masse:

40 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

RASSF7, also named as HRAS1 and HRC1, is a member of the RASSF family contain a highly conserved Ras association domain (RalGDS/AF-6 or RA) and function as Ras effectors tumor suppressors. RASSF7, RASSF8, RASSF9 and RASSF10 are four evolutionary conserved proteins that are collectively known as the N-terminal RASSF family. RASSF7 is believed to play an essential role in mitosis. RASSF7 has been shown to influence cell failure to form mitotic spindle and also cause nuclear breakdown and apoptosis. The antibody can recognizes all the 3 isoform of RASSF7.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

