

Nur für Forschungszwecke

CYP4F12-Specific Polyklonaler Antikörper

Katalog-Nr.: 20013-1-AP



Allgemeine Informationen

Katalog-Nr.:

20013-1-AP

Größe:

150ul, Konzentration: 253 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_023944

GeneID (NCBI):

66002

Vollständiger Name:

cytochrome P450, family 4, subfamily
F, polypeptide 12

Berechnete Masse:

60 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

CYP4F12 belongs to the cytochrome P450 family. Cytochromes P450 are a group of heme-thiolate monooxygenases. CYP4F12 catalyzes leukotriene B₄ omega-hydroxylation and arachidonic acid omega-hydroxylation but with an activity much lower than that of CYP4F2. It catalyzes the hydroxylation of the antihistamine ebastine. This antibody is specific to CYP4F12. It has no cross reaction to other CYP4F member.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

