

Nur für Forschungszwecke

# MED13-Specific Polyklonaler Antikörper

Katalog-Nr.: 20029-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
20029-1-AP

Größe:  
150ul, Konzentration: 293 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

GenBank-Zugangsnummer:  
NM\_005121

GeneID (NCBI):  
9969

Vollständiger Name:  
mediator complex subunit 13

Berechnete Masse:  
239 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human, Maus, Ratte

## Hintergrundinformationen

MED13, also named as ARC250, KIAA0593, THRAP1 and TRAP240, belongs to the Mediator complex subunit 13 family. MED13 is a component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional preinitiation complex with RNA polymerase II and the general transcription factors. This antibody is specific to MED13. It has no cross reaction to MED13L.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

