

Nur für Forschungszwecke

# APOL5-Specific Polyklonaler Antikörper

Katalog-Nr.: 20068-1-AP



## Allgemeine Informationen

Katalog-Nr.:

20068-1-AP

Größe:

150ul, Konzentration: 260 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM\_030642

GeneID (NCBI):

80831

Vollständiger Name:

apolipoprotein L, 5

Berechnete Masse:

47 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

## Hintergrundinformationen

APOL5 belongs to the apolipoprotein L family. The pathological expression of the apolipoproteins during pregnancy is implicated in fetal growth retardation, preeclampsia, and the onset of adult atherosclerosis [PubMed: 11374903]. APOL5 may affect the movement of lipids in the cytoplasm or allow the binding of lipids to organelles. This antibody is specific to APOL5.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

