

Nur für Forschungszwecke

GLS2 Polyklonaler Antikörper

Katalog-Nr.: 20171-1-AP



Allgemeine Informationen

Katalog-Nr.: 20171-1-AP	GenBank-Zugangsnummer: BC059370	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 333 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 27165	
Wirt: Kaninchen	Vollständiger Name: glutaminase 2 (liver, mitochondrial)	
Isotyp: IgG	Berechnete Masse: 602 aa, 66 kDa	
Immunogen Katalognummer: AG14003		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

GLS2, also named as GA and GLS, belongs to the glutaminase family. It is Glutaminase liver isoform . GLS2 plays an important role in the regulation of glutamine catabolism. It promotes mitochondrial respiration and increases ATP generation in cells by catalyzing the synthesis of glutamate and alpha-ketoglutarate. GLS2 may play a role in preventing tumor proliferation. This antibody is specific to GLS2.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

