

Nur für Forschungszwecke

# SRGAP1-Specific Polyklonaler Antikörper

Katalog-Nr.: 20222-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
20222-1-AP

Größe:  
150ul, Konzentration: 280 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

GenBank-Zugangsnummer:  
NM\_020762

GeneID (NCBI):  
57522

Vollständiger Name:  
SLIT-ROBO Rho GTPase activating  
protein 1

Berechnete Masse:  
124 kDa

Reinigungsmethode:  
Antigen-Affinitätsreinigung

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human, Maus, Ratte

## Hintergrundinformationen

SRGAP1, also named as ARHGAP13 and KIAA1304, is a GTPase-activating protein for RhoA and Cdc42 small GTPases. Together with CDC42, SRGAP1 seems to be involved in the pathway mediating the repulsive signaling of Robo and Slit proteins in neuronal migration. SLIT2, probably through interaction with ROBO1, increases the interaction of SRGAP1 with ROBO1 and inactivates CDC42. This antibody is specific to SRGAP1. It can't recognize SRGAP2 and SRGAP3.

## Lagerung

Lagerungsbedingungen:  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
Lagerungspuffer:  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

