

Nur für Forschungszwecke

TSPAN4 Polyklonaler Antikörper

Katalog-Nr.: 20318-1-AP



Allgemeine Informationen

Katalog-Nr.: 20318-1-AP	GenBank-Zugangsnummer: BC000389	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 213 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 7106	
Wirt: Kaninchen	Vollständiger Name: tetraspanin 4	
Isotyp: IgG	Berechnete Masse: 238 aa, 26 kDa	
Immunogen Katalognummer: AG14172		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

TSPAN4 (also known as NAG2 or TM4SF7) is a membrane protein of the tetraspanin superfamily, which are characterized by the presence of four conserved transmembrane regions. Tetraspanins have been involved in diverse processes such as cell activation and proliferation, adhesion and motility, differentiation, and cancer. TSPAN4 forms complexes with integrins and other tetraspanins (PMID: 9360996).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

