

Nur für Forschungszwecke

ABCB5 Polyklonaler Antikörper

Katalog-Nr.:20703-1-AP



Allgemeine Informationen

Katalog-Nr.:
20703-1-AP

Größe:
150ul , Konzentration: 307 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

GenBank-Zugangsnummer:
NM_178559

GeneID (NCBI):
340273

Vollständiger Name:
ATP-binding cassette, sub-family B
(MDR/TAP), member 5

Berechnete Masse:
90 kDa

Beobachtete Masse:
89-110 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB : A431-Zellen, humanes Hodengewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



A431 cells were subjected to SDS PAGE followed by western blot with 20703-1-AP (ABCB5 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.