

Nur für Forschungszwecke

ADAM29 Polyklonaler Antikörper

Katalog-Nr.: 20736-1-AP



Allgemeine Informationen

Katalog-Nr.:

20736-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

GenBank-Zugangsnummer:

NM_014269

GeneID (NCBI):

11086

Vollständiger Name:

ADAM metallopeptidase domain 29

Berechnete Masse:

93 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

ADAM29 is also named as CT73. It may be involved in spermatogenesis and fertilization. ADAM29 is processed to mature forms (70-85 kDa) between the stages of testicular sperm and cauda epididymal sperm. The cytoplasmic domain of ADAM29 can be processed or modified in the mature form so that the antibody might be unable to recognize the protein (PMID:19129510). It has 3 isoforms produced by alternative splicing. The full length protein has a signal peptide, a propeptide and 9 glycosylation sites. This antibody is specific to ADAM29.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

