

## Allgemeine Informationen

<b>Katalog-Nr.:</b> 20875-1-AP	<b>GenBank-Zugangsnummer:</b> BC011698	<b>Reinigungsmethode:</b> Antigen-Affinitätsreinigung
<b>Größe:</b> 150ul, Konzentration: 550 µg/ml von Nanodrop und 387 µg/ml durch die Bradford-Methode mit BSA als Standard;	<b>GeneID (NCBI):</b> 54902	<b>Empfohlene Verdünnungen:</b> WB 1:500-1:1000 IHC 1:50-1:500
<b>Wirt:</b> Kaninchen	<b>Vollständiger Name:</b> tetra(tryptophan) repeat domain 19	
<b>Isotyp:</b> IgG	<b>Berechnete Masse:</b> 501 aa, 56 kDa	
<b>Immunogen Katalognummer:</b> AG14974	<b>Beobachtete Masse:</b> 32-35 kDa	

## Anwendungen

<b>Geprüfte Anwendungen:</b> IHC, WB, ELISA	<b>Positivkontrollen:</b>
<b>Getestete Reaktivität:</b> Human, Maus	<b>WB:</b> Mausherzgewebe, HeLa-Zellen, Rattenherzgewebe
<b>Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.</b>	<b>IHC:</b> humanes Mammakarzinomgewebe, Mausherzgewebe

## Hintergrundinformationen


TTC19 is a subunit of mitochondrial respiratory chain complex III and is proposed to preserve the structural and functional integrity of complex III. Loss-of-function mutations in TTC19 have been found in patients with complex III deficiency associated with heterogeneous neurological syndromes. The precursor TTC19 was detected as a 43 kDa cytoplasmic protein, while the mature TTC19 with the mitochondrial targeting sequence cleaved was observed as 35 kDa by SDS-PAGE. (PMID: 21278747)

## Lagerung


**Lagerungsbedingungen:**  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
**Lagerungspuffer:**  
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.  
**Aliquotieren ist nicht notwendig bei -20°C Lagerung**

\*\*\* 20ul-Größen enthalten 0.1% BSA

## Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 20875-1-AP (TTC19 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue slide using 20875-1-AP (TTC19 antibody) at dilution of 1:200 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0)).