

Nur für Forschungszwecke

# ASTL Polyklonaler Antikörper

Katalog-Nr.: 21418-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
21418-1-AP

Größe:

150ul, Konzentration: 240 µg/ml von  
Nanodrop und 200 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:  
AG15997

GenBank-Zugangsnummer:  
BC107127

GeneID (NCBI):  
431705

Vollständiger Name:  
astacin-like metallo-endopeptidase  
(M12 family)

Berechnete Masse:  
431 aa, 46 kDa

Beobachtete Masse:  
43 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:  
WB 1:200-1:1000

## Anwendungen

Geprüfte Anwendungen:  
WB, ELISA

Getestete Reaktivität:  
Human

Positivkontrollen:

WB: HEK-293-Zellen,

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



HEK-293 cells were subjected to SDS PAGE followed by western blot with 21418-1-AP (ASTL Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.