

Nur für Forschungszwecke

ESYT1 Polyklonaler Antikörper

Katalog-Nr.: 21450-1-AP



Allgemeine Informationen

Katalog-Nr.: 21450-1-AP	GenBank-Zugangsnummer: BC004998	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 300 µg/ml von Nanodrop;	GeneID (NCBI): 23344	Empfohlene Verdünnungen: WB 1:2000-1:8000
Wirt: Kaninchen	Vollständiger Name: family with sequence similarity 62 (C2 domain containing), member A	
Isotyp: IgG	Berechnete Masse: 1104 aa, 123 kDa	
Immunogen Katalognummer: AG13895	Beobachtete Masse: 123 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : HEK-293T-Zellen, HeLa-Zellen, Jurkat-Zellen
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 21450-1-AP (ESYT1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.