

Nur für Forschungszwecke

TTC31 Polyklonaler Antikörper

Katalog-Nr.: 21693-1-AP



Allgemeine Informationen

| | | |
|--|---|---|
| Katalog-Nr.: 21693-1-AP | GenBank-Zugangsnummer: BC098365 | Reinigungsmethode: Antigen-Affinitätsreinigung |
| Größe: 150ul, Konzentration: 450 µg/ml von Nanodrop und 227 µg/ml durch die Bradford-Methode mit BSA als Standard; | GeneID (NCBI): 64427 | Empfohlene Verdünnungen: WB 1:500-1:2000 IHC 1:20-1:200 IF 1:20-1:200 |
| Wirt: Kaninchen | Vollständiger Name: tetra(Tricopeptide repeat domain 31 | |
| Isotyp: IgG | Berechnete Masse: 358 aa, 40 kDa | |
| Immunogen Katalognummer: AG16281 | Beobachtete Masse: 57 kDa | |

Anwendungen

| | |
|---|--|
| Geprüfte Anwendungen: IF, IHC, WB, ELISA | Positivkontrollen: WB: PC-3-Zellen, IHC: humanes Hodengewebe, humanes Lebergewebe IF: PC-3-Zellen, |
| Getestete Reaktivität: Human | |
| Hinweis-IHC: Antigenmaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen. | |

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




PC-3 cells were subjected to SDS PAGE followed by western blot with 21693-1-AP (TTC31 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human testis using 21693-1-AP (TTC31 antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis using 21693-1-AP (TTC31 antibody) at dilution of 1:50 (under 40x lens).



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using 21693-1-AP (TTC31 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).