

Nur für Forschungszwecke

BMT2 Polyklonaler Antikörper

Katalog-Nr.: 21744-1-AP



Allgemeine Informationen

Katalog-Nr.: 21744-1-AP	GenBank-Zugangsnummer: BC114615	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 300 µg/ml von Nanodrop;	GeneID (NCBI): 154743	Empfohlene Verdünnungen: WB 1:1000-1:6000
Wirt: Kaninchen	Vollständiger Name: chromosome 7 open reading frame 60	
Isotyp: IgG	Berechnete Masse: 405 aa, 46 kDa	
Immunogen Katalognummer: AG16476	Beobachtete Masse: 46-50 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : HeLa-Zellen, U-87 MG-Zellen
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 21744-1-AP (C7orf60 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.