

Nur für Forschungszwecke

KIAA0895 Polyklonaler Antikörper

Katalog-Nr.: 21930-1-AP



Allgemeine Informationen

Katalog-Nr.:
21930-1-AP

Größe:
150ul, Konzentration: 700 µg/ml von
Nanodrop und 353 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG16533

GenBank-Zugangsnummer:
BC110422

GeneID (NCBI):
23366
Vollständiger Name:
KIAA0895

Berechnete Masse:
520 aa, 61 kDa

Beobachtete Masse:
45 kDa, 60 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:

WB: Maushodengewebe, HEK-293-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse testis tissue were subjected to SDS PAGE followed by western blot with 21930-1-AP (KIAA0895 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.