

Nur für Forschungszwecke

SOX2 Polyklonaler Antikörper

Katalog-Nr.:21965-1-AP



Allgemeine Informationen

Katalog-Nr.: 21965-1-AP	GenBank-Zugangsnummer: BC013923	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul	GeneID (NCBI): 6657	
Wirt: Kaninchen	Vollständiger Name: SRY (sex determining region Y)-box 2	
Isotyp: IgG	Berechnete Masse: 317 aa, 34 kDa	
Immunogen Katalognummer: AG16794		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Sox2, also known as SRY (sex determining region Y)-box 2, is a transcription factor essential for maintaining self-renewal of undifferentiated ES cells and is one of the key transcription factors used to reprogram mouse and human fibroblasts to a pluripotent state. Sox2 expressed in undifferentiated pluripotent stem cells and germ cells during development. Affinity purified rabbit anti-Sox2 antibody can be used to demonstrate pluripotency of ES and iPS cells. This antibody is a rabbit polyclonal antibody raised against residues near the C terminus of human SOX2.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

