

Nur für Forschungszwecke

PPARD Polyklonaler Antikörper

Katalog-Nr.: 21986-1-AP



Allgemeine Informationen

Katalog-Nr.: 21986-1-AP	GenBank-Zugangsnummer: BC002715	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul	GeneID (NCBI): 5467	
Wirt: Kaninchen	Vollständiger Name: peroxisome proliferator-activated receptor delta	
Isotyp: IgG	Berechnete Masse: 441 aa, 50 kDa	
Immunogen Katalognummer: AG17135		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Peroxisome proliferator-activator receptor delta (PPARD), a ligand-activated transcription factor, belongs to the PPARs family, which get its name for their chemicals that induce proliferation of peroxisomes, organelles that contribute to the oxidation of fatty acids. PPARD is a receptor for peroxisome proliferator and once activated by a ligand, it binds to and activates the promoter elements of target genes, such as acylCoA oxidase gene. But it act as a repressor for the NPC1L1 gene.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

