

DDX56 Polyklonaler Antikörper

Katalog-Nr.: 21989-1-AP

Allgemeine Informationen

Katalog-Nr.:

21989-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von Nanodrop und 413 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG17272

GenBank-Zugangsnummer:

BC001235

GeneID (NCBI):

54606

Vollständiger Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 56

Berechnete Masse:

62 kDa

Beobachtete Masse:

62 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : HepG2-Zellen, Maushirngewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HepG2 cells were subjected to SDS PAGE followed by western blot with 21989-1-AP (DDX56 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.