

Nur für Forschungszwecke

Cadherin-20 Polyklonaler Antikörper

Katalog-Nr.: 22049-1-AP



Allgemeine Informationen

Katalog-Nr.: 22049-1-AP	GenBank-Zugangsnummer: BC101048	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 750 µg/ml von Nanodrop und 440 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 28316	Empfohlene Verdünnungen: WB 1:500-1:2000 IHC 1:2000-1:8000
Wirt: Kaninchen	Vollständiger Name: cadherin 20, type 2	
Isotyp: IgG	Berechnete Masse: 801 aa, 89 kDa	
Immunogen Katalognummer: AG16990	Beobachtete Masse: 120-130 kDa	

Anwendungen

Geprüfte Anwendungen: IHC, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human, Maus, Ratte	WB: humanes Plazenta-Gewebe, HeLa-Zellen, Maushirngewebe
Hinweis-IHC: Antigenmaskierung mit TE- Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigenmaskierung auch mit Citratpuffer pH 6,0 erfolgen.	IHC: humanes Kolongewebe, Maushirngewebe, Rattenhirngewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



Human placenta tissue were subjected to SDS PAGE followed by western blot with 22049-1-AP (Cadherin-20 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 22049-1-AP (Cadherin-20 antibody) at dilution of 1:4000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon tissue slide using 22049-1-AP (Cadherin-20 antibody) at dilution of 1:4000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).