

Nur für Forschungszwecke

MUSK Polyklonaler Antikörper

Katalog-Nr.: 22112-1-AP



Allgemeine Informationen

Katalog-Nr.:
22112-1-AP

Größe:
150ul, Konzentration: 453 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG17371

GenBank-Zugangsnummer:
BC109098

GeneID (NCBI):
4593

Vollständiger Name:
muscle, skeletal, receptor tyrosine
kinase

Berechnete Masse:
869 aa, 97 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB: Jurkat-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Jurkat cells were subjected to SDS PAGE followed by western blot with 22112-1-AP (MUSK Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.