

HOXC11 Polyklonaler Antikörper

Katalog-Nr.: 22118-1-AP

Allgemeine Informationen

Katalog-Nr.:
22118-1-AP

Größe:

150ul, Konzentration: 500 µg/ml von
Nanodrop und 307 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG17806

GenBank-Zugangsnummer:

BC001543

GeneID (NCBI):

3227

Vollständiger Name:

homeobox C11

Berechnete Masse:

304 aa, 34 kDa

Beobachtete Masse:

47 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

HOXC11, also named as Homeobox protein Hox 3H, is a 304 amino acid protein, which contains one homeobox DNA-binding domain and belongs to the Abd-B homeobox family. HOXC11 is a sequence specific transcription factor, which is a part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

