

Nur für Forschungszwecke

CACNA2D2 Polyklonaler Antikörper

Katalog-Nr.: 22145-1-AP



Allgemeine Informationen

Katalog-Nr.:

22145-1-AP

Größe:

150ul

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG17580

GenBank-Zugangsnummer:

BC152438

GeneID (NCBI):

9254

Vollständiger Name:

calcium channel, voltage-dependent,
alpha 2/delta subunit 2

Berechnete Masse:

1145 aa, 129 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Maushirngewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse brain tissue were subjected to SDS PAGE
followed by western blot with 22145-1-AP
(CACNA2D2 Antibody) at dilution of 1:600
incubated at room temperature for 1.5 hours.