

Nur für Forschungszwecke

SOSTDC1 Polyklonaler Antikörper

Katalog-Nr.:22213-1-AP



Allgemeine Informationen

Katalog-Nr.:
22213-1-AP

Größe:
150ul , Konzentration: 280 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG17538

GenBank-Zugangsnummer:
BC008484

GeneID (NCBI):
25928

Vollständiger Name:
sclerostin domain containing 1

Berechnete Masse:
206 aa, 23 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human, Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 22213-1-AP (SOSTDC1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.