

Nur für Forschungszwecke

HOXB4 Polyklonaler Antikörper

Katalog-Nr.:22256-1-AP



Allgemeine Informationen

Katalog-Nr.:
22256-1-AP

Größe:
150ul, Konzentration: 650 µg/ml von
Nanodrop und 280 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG17761

GenBank-Zugangsnummer:
BC049204

GeneID (NCBI):
3214

Vollständiger Name:
homeobox B4

Berechnete Masse:
251 aa, 28 kDa

Beobachtete Masse:
28-35 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB: Jurkat-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Jurkat cells were subjected to SDS PAGE followed by western blot with 22256-1-AP (HOXB4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.