

Nur für Forschungszwecke

CDX4 Polyklonaler Antikörper

Katalog-Nr.:22943-1-AP



Allgemeine Informationen

Katalog-Nr.:
22943-1-AP

Größe:

150ul, Konzentration: 450 µg/ml von
Nanodrop und 260 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG19078

GenBank-Zugangsnummer:

BC128233

GeneID (NCBI):

1046

Vollständiger Name:

caudal type homeobox 4

Berechnete Masse:

284 aa, 30 kDa

Beobachtete Masse:

37 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : SKOV-3-Zellen, Maus-Thymusgewebe, Sp2/O-
Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



SKOV-3 cells were subjected to SDS PAGE followed by western blot with 22943-1-AP (CDX4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.