

Nur für Forschungszwecke

DSG4 Polyklonaler Antikörper

Katalog-Nr.: 23087-1-AP



Allgemeine Informationen

Katalog-Nr.:

23087-1-AP

Größe:

150ul, Konzentration: 1000 µg/ml von 147409

Nanodrop und 633 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG19520

GenBank-Zugangsnummer:

BC132907

GeneID (NCBI):

147409

Vollständiger Name:

desmoglein 4

Berechnete Masse:

1040 aa, 114 kDa

Beobachtete Masse:

160 kDa

Reinigungsmethode:

Antigen-affinitätsgereinigt

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Mauhautgewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse skin tissue were subjected to SDS PAGE followed by western blot with 23087-1-AP (DSG4 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.