

Nur für Forschungszwecke

HDAC9 Polyklonaler Antikörper

Katalog-Nr.: 23328-1-AP



Allgemeine Informationen

Katalog-Nr.:
23328-1-AP

Größe:
150ul, Konzentration: 433 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG19930

GenBank-Zugangsnummer:
BC152405

GeneID (NCBI):
9734

Vollständiger Name:
histone deacetylase 9

Berechnete Masse:
1011 aa, 111 kDa

Beobachtete Masse:
100 kDa

Reinigungsmethode:
Antigen-affinitätsgereinigt

Anwendungen

Geprüfte Anwendungen:
IHC, ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 23328-1-AP (HDAC9 Antibody) at dilution of 1:200 (under 40x lens).