

Nur für Forschungszwecke

SCGF Polyklonaler Antikörper

Katalog-Nr.: 24138-1-AP



Allgemeine Informationen

Katalog-Nr.: 24138-1-AP	GenBank-Zugangsnummer: BC005810	Reinigungsmethode: Antigen-affinitätsgereinigt
Größe: 150ul, Konzentration: 340 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 6320	
Wirt: Kaninchen	Vollständiger Name: C-type lectin domain family 11, member A	
Isotyp: IgG	Berechnete Masse: 36 kDa	
Immunogen Katalognummer: AG21365		

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

CLEC11A, also named as CLECSF3, LSLCL, SCGF and p47, stimulates the proliferation and differentiation of hematopoietic precursor cells from various lineages, including erythrocytes, lymphocytes, granulocytes and macrophages. CLEC11A acts synergistically with other cytokines, including IL-3, GCSF, GMCSF and FLT3 ligand. It suppresses SCF-stimulated erythrocyte proliferation.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

