

Nur für Forschungszwecke

# Shootin 1 Polyklonaler Antikörper

Katalog-Nr.: 24383-1-AP

Vorgestelltes Produkt



## Allgemeine Informationen

Katalog-Nr.:  
24383-1-AP

Größe:  
150ul, Konzentration: 350 µg/ml von  
Nanodrop und 167 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG18294

GenBank-Zugangsnummer:  
BC022348

GeneID (NCBI):  
57698

Vollständiger Name:  
KIAA1598

Berechnete Masse:  
72 kDa

Beobachtete Masse:  
60-75 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

## Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: Neuro-2a-Zellen, Rattenhirngewebe

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



WB result of Shootin 1 antibody (24383-1-AP; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-Shootin 1 transfected Neuro-2a cells.