

Allgemeine Informationen

Katalog-Nr.:

24426-1-AP

Größe:

150ul, Konzentration: 400 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG18833

GenBank-Zugangsnummer:

BC146464

GeneID (NCBI):

128859

Vollständiger Name:

bactericidal/permeability-increasing
protein-like 3

Berechnete Masse:

453 aa, 50 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

BPIL3 (Bactericidal/permeability-increasing protein-like 3), also known as BPIFB6 (BPI fold-containing family B member 6), belongs to the lipid transfer (LT)/lipopolysaccharide (LPS)-binding protein family. LT/LBP proteins are capable of binding phospholipids and LPS. Expression of BPIL3 gene was detected at very low levels in normal tonsils, and at high levels in hypertrophic tonsils (PMID: 12185532). BPIL3 may play a role in innate immunity.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

