

Nur für Forschungszwecke

ZNF284 Polyklonaler Antikörper

Katalog-Nr.:24681-1-AP



Allgemeine Informationen

Katalog-Nr.:
24681-1-AP

Größe:

150ul, Konzentration: 350 µg/ml von
Nanodrop und 193 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:
AG20336

GenBank-Zugangsnummer:

BC156663

GeneID (NCBI):

342909

Vollständiger Name:

zinc finger protein 284

Berechnete Masse:

593 aa, 69 kDa

Beobachtete Masse:

69 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB: Mausebergewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



mouse liver tissue were subjected to SDS PAGE followed by western blot with 24681-1-AP (ZNF284 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.