

Nur für Forschungszwecke

ZNF674 Polyklonaler Antikörper

Katalog-Nr.: 24747-1-AP



Allgemeine Informationen

Katalog-Nr.:

24747-1-AP

Größe:

150ul, Konzentration: 500 µg/ml von Nanodrop und 300 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG21576

GenBank-Zugangsnummer:

BC157038

GeneID (NCBI):

641339

Vollständiger Name:

zinc finger family member 674

Berechnete Masse:

581 aa, 67 kDa

Beobachtete Masse:

63-67 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : COLO 320-Zellen, HepG2-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 24747-1-AP (ZNF674 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.