

Nur für Forschungszwecke

FOXD4 Polyklonaler Antikörper

Katalog-Nr.: 24835-1-AP



Allgemeine Informationen

Katalog-Nr.:
24835-1-AP

Größe:
150ul, Konzentration: 1000 µg/ml von 2298
Nanodrop und 440 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG19226

GenBank-Zugangsnummer:
BC136570

GeneID (NCBI):

Vollständiger Name:
forkhead box D4

Berechnete Masse:
439 aa, 47 kDa

Beobachtete Masse:
65-70 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:2000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

WB: NIH/3T3-Zellen, Jurkat-Zellen, RAW 264.7-Zellen,
U-937-Zellen

Hintergrundinformationen

Forkhead box protein D4 (FOXD4) also named as FKHL9 or Myeloid factor alpha is a 439 amino acid protein, which contains one fork head DNA-binding domain. FOXD4L6 localizes in the nucleus. FOXD4 as an embryonic transcriptional regulator involves in leukemogenesis. The molecular weight of FOXD4L6 is 47kDa, but our experiment always detected a 60-70kDa protein by western blot.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 24835-1-AP (FOXD4 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.