

Nur für Forschungszwecke

POTEH Polyklonaler Antikörper

Katalog-Nr.:24837-1-AP



Allgemeine Informationen

Katalog-Nr.: 24837-1-AP	GenBank-Zugangsnummer: BC148758	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , Konzentration: 280 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 23784	Empfohlene Verdünnungen: WB 1:200-1:1000
Wirt: Kaninchen	Vollständiger Name: POTE ankyrin domain family, member H	
Isotyp: IgG	Berechnete Masse: 545 aa, 61 kDa	
Immunogen Katalognummer: AG20275	Beobachtete Masse: 61 kDa	

Anwendungen

Geprüfte Anwendungen: WB,ELISA	Positivkontrollen: WB : DU 145-Zellen,
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



DU 145 cells were subjected to SDS PAGE followed by western blot with 24837-1-AP (POTEH Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.