

Nur für Forschungszwecke

# CNTD1 Polyklonaler Antikörper

Katalog-Nr.:24854-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
24854-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von  
Nanodrop und 247 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:  
AG21338

GenBank-Zugangsnummer:

BC026187

GeneID (NCBI):

124817

Vollständiger Name:

cyclin N-terminal domain containing  
1

Berechnete Masse:

330 aa, 37 kDa

Beobachtete Masse:

35-37 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

## Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : Maus-Kolongewebe,

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



mouse colon tissue were subjected to SDS PAGE followed by western blot with 24854-1-AP (CNTD1 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.