

Nur für Forschungszwecke

ANKFY1 Polyklonaler Antikörper

Katalog-Nr.: 24890-1-AP

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:

24890-1-AP

Größe:

150ul, Konzentration: 320 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG19267

GenBank-Zugangsnummer:

BC152991

GeneID (NCBI):

51479

Vollständiger Name:

ankyrin repeat and FYVE domain
containing 1

Berechnete Masse:

1169 aa, 129 kDa

Beobachtete Masse:

105 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus

Hintergrundinformationen

ANKFY1 also named as ANKHZN or Ankyrin repeats hooked to a zinc finger motif is a 1169 amino acid protein, which is a cytoplasm protein on embryonic stem cells. ANKFY1 protein is ubiquitously expressed in a spatiotemporal specific manner and is located on endosomes. ANKFY1 is expressed at various levels in human tissues and also present in both membrane and soluble fractions obtained on subcellular fractionation. ANKFY1 is a single copy gene consisting of predicted 25 exons in the human genome, and has been mapped to human chromosome 17p13.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Chao Weng	28123425	Neural Regen Res	

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

