

Nur für Forschungszwecke

# BTNL8 Polyklonaler Antikörper

Katalog-Nr.:24959-1-AP



## Allgemeine Informationen

Katalog-Nr.:  
24959-1-AP

Größe:  
150ul, Konzentration: 700 µg/ml von  
Nanodrop und 347 µg/ml durch die  
Bradford-Methode mit BSA als  
Standard;

Wirt:  
Kaninchen

Isotyp:  
IgG

Immunogen Katalognummer:  
AG21203

GenBank-Zugangsnummer:  
BC119696

GeneID (NCBI):  
79908

Vollständiger Name:  
butyrophilin-like 8

Berechnete Masse:  
500 aa, 57 kDa

Beobachtete Masse:  
50 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

## Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: MCF-7-Zellen,

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



MCF-7 cells were subjected to SDS PAGE followed by western blot with 24959-1-AP (BTNL8 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.