

SP2 Polyklonaler Antikörper

Katalog-Nr.: 25000-1-AP

Allgemeine Informationen

Katalog-Nr.:
25000-1-AP

Größe:
150ul, Konzentration: 500 µg/ml von
Nanodrop und 300 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG8060

GenBank-Zugangsnummer:
BC016680

GeneID (NCBI):
6668

Vollständiger Name:
Sp2 transcription factor

Berechnete Masse:
65 kDa

Beobachtete Masse:
72 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:
WB, ELISA

Getestete Reaktivität:
Human, Ratte

Positivkontrollen:

WB: Jurkat-Zellen, C6-Zellen, HEK-293-Zellen, K-562-Zellen

Hintergrundinformationen

SP2 also named as Sp2 transcription factor or KIAA0048 is a 613 amino acid protein, which contains 3 C2H2-type zinc fingers and belongs to the sp1 C2H2-type zinc finger protein family. SP2 binds to box promoters elements and selectively activates mRNA synthesis from genes that contain functional recognition sites in nucleus. SP2 may regulate expression of the T-cell antigen receptor (TCR) gene.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Jurkat cells were subjected to SDS PAGE followed by western blot with 25000-1-AP (SP2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



K-562 cells were subjected to SDS PAGE followed by western blot with 25000-1-AP (SP2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.