

Nur für Forschungszwecke

POU3F4 Polyklonaler Antikörper

Katalog-Nr.:25114-1-AP



Allgemeine Informationen

Katalog-Nr.:
25114-1-AP

Größe:
150ul , Konzentration: 500 µg/ml von
Nanodrop und 220 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG18925

GenBank-Zugangsnummer:
BC146551

GeneID (NCBI):
5456

Vollständiger Name:
POU class 3 homeobox 4

Berechnete Masse:
361 aa, 39 kDa

Beobachtete Masse:
40-45 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:200-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:

WB : fetales humanes Hirngewebe,

Hintergrundinformationen

OTF9 also named as BRN4 and POU3F4, belongs to the POU transcription factor family and Class-3 subfamily. OTF9 is a transcription factor which exerts its primary action widely during early neural development and in a very limited set of neurons in the mature brain. Mutation of OTF9 will cause of X-linked deafness type 3 (DFN3).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 25114-1-AP (POU3F4 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.