

Allgemeine Informationen

Katalog-Nr.:

25216-1-AP

Größe:

150ul, Konzentration: 400 µg/ml von Nanodrop und 260 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG17165

GenBank-Zugangsnummer:

BC099661

GeneID (NCBI):

7678

Vollständiger Name:

zinc finger protein 124

Berechnete Masse:

351 aa, 40 kDa

Beobachtete Masse:

40-45 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:1000-1:4000

IP 0.5-4.0 µg für IP und 1:1000-1:4000

für WB

IHC 1:20-1:200

IF 1:10-1:100

Anwendungen

Geprüfte Anwendungen:

IF, IHC, IP, WB, ELISA

Getestete Reaktivität:

Human

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : SH-SY5Y-Zellen, U-937-Zellen

IP : SH-SY5Y-Zellen,

IHC : humanes Gliomgewebe,

IF : HepG2-Zellen,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 25216-1-AP (ZNF 124 Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



IP Result of anti-ZNF124 (IP:25216-1-AP, 4ug; Detection:25216-1-AP 1:2000) with SH-SY5Y cells lysate 1600ug.



Immunohistochemical analysis of paraffin-embedded human gliomas tissue slide using 25216-1-AP (ZNF124 Antibody) at dilution of 1:50 (under 40x lens).



Immunofluorescent analysis of HepG2 cells using 25216-1-AP (ZNF124 antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).