

Nur für Forschungszwecke

PDZD8 Polyklonaler Antikörper

Katalog-Nr.: 25512-1-AP



Allgemeine Informationen

Katalog-Nr.: 25512-1-AP	GenBank-Zugangsnummer: BC028375	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul, Konzentration: 700 µg/ml von Nanodrop;	GeneID (NCBI): 118987	Empfohlene Verdünnungen: WB 1:500-1:1000
Wirt: Kaninchen	Vollständiger Name: PDZ domain containing 8	
Isotyp: IgG	Berechnete Masse: 1154 aa, 129 kDa	
Immunogen Katalognummer: AG21914	Beobachtete Masse: 150-160 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : HCT 116-Zellen, HEK-293-Zellen, HeLa-Zellen
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 25512-1-AP (PDZD8 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.