

Nur für Forschungszwecke

ZC3H7B Polyklonaler Antikörper

Katalog-Nr.:25624-1-AP



Allgemeine Informationen

Katalog-Nr.:
25624-1-AP

Größe:
150ul, Konzentration: 650 µg/ml von
Nanodrop und 327 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG22555

GenBank-Zugangsnummer:
BC152558

GeneID (NCBI):
23264

Vollständiger Name:
zinc finger CCCH-type containing 7B

Berechnete Masse:
993 aa, 112 kDa

Beobachtete Masse:
111 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human, Maus

Positivkontrollen:

WB : HUVEC-Zellen, Maushirngewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



HUVEC cells were subjected to SDS PAGE followed by western blot with 25624-1-AP (ZC3H7B Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.