

Nur für Forschungszwecke

ZNF668 Polyklonaler Antikörper

Katalog-Nr.:25669-1-AP



Allgemeine Informationen

Katalog-Nr.:
25669-1-AP

Größe:
150ul, Konzentration: 400 µg/ml von
Nanodrop und 273 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG22506

GenBank-Zugangsnummer:
BC021997

GeneID (NCBI):
79759

Vollständiger Name:
zinc finger protein 668

Berechnete Masse:
619 aa, 68 kDa

Beobachtete Masse:
70 kDa

Reinigungsmethode:
Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:
WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human

Positivkontrollen:

WB : COLO 320-Zellen, A549-Zellen, HepG2-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



COLO 320 cells were subjected to SDS PAGE followed by western blot with 25669-1-AP (ZNF668 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.