

PLA2R1 Polyklonaler Antikörper

Katalog-Nr.: 25695-1-AP

Allgemeine Informationen

Katalog-Nr.:

25695-1-AP

Größe:

150ul, Konzentration: 213 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG22448

GenBank-Zugangsnummer:

BC140823

GeneID (NCBI):

22925

Vollständiger Name:

phospholipase A2 receptor 1, 180kDa

Berechnete Masse:

1463 aa, 169 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Secretory phospholipase A2 (sPLA2) are implicated in a number of crucial physiological and pathological processes in mammals. Different types of membrane receptors for sPLA2s have been identified. The M-type receptor, PLA2R1 (phospholipase A2 receptor 1, 180 kDa), is a type I transmembrane glycoprotein containing extracellularly a cysteine-rich domain, a fibronectin type II domain and eight carbohydrate-recognition domains. It is a receptor for phospholipase sPLA2-IB/PLA2G1B but not sPLA2-IIA/PLA2G2A. It is also able to bind to snake PA2-like toxins. PLA2R1 may be involved in responses in proinflammatory cytokine productions during endotoxic shock. It also has endocytic properties and rapidly internalizes sPLA2 ligands, which is particularly important for the clearance of extracellular sPLA2s to protect their potent enzymatic activities. It has been recently shown that PLA2R1 is the major target antigen in idiopathic membranous nephropathy (iMN). Patients with iMN have antibodies against a conformation-dependent epitope in PLA2R1.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

