

Nur für Forschungszwecke

ZNF419 Polyklonaler Antikörper

Katalog-Nr.: 25724-1-AP



Allgemeine Informationen

Katalog-Nr.:

25724-1-AP

Größe:

150ul, Konzentration: 650 µg/ml von Nanodrop und 340 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG22802

GenBank-Zugangsnummer:

BC008297

GeneID (NCBI):

79744

Vollständiger Name:

zinc finger protein 419

Berechnete Masse:

510 aa, 59 kDa

Beobachtete Masse:

59 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB: A549-Zellen, Jurkat-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



A549 cells were subjected to SDS PAGE followed by western blot with 25724-1-AP (ZNF419 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.