

Nur für Forschungszwecke

ZMYM3 Polyklonaler Antikörper

Katalog-Nr.:25742-1-AP



Allgemeine Informationen

Katalog-Nr.:
25742-1-AP

Größe:
150ul , Konzentration: 500 µg/ml von
Nanodrop und 333 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:
Kaninchen

Isotyp:
IgG

Immunogen Katalognummer:
AG22701

GenBank-Zugangsnummer:
BC013009

GeneID (NCBI):
9203

Vollständiger Name:
zinc finger, MYM-type 3

Berechnete Masse:
1370 aa, 152 kDa

Beobachtete Masse:
190-200 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Empfohlene Verdünnungen:

WB 1:500-1:1000
IP 0.5-4.0 ug für IP und 1:200-1:1000
für WB

Anwendungen

Geprüfte Anwendungen:

IP, WB,ELISA

Getestete Reaktivität:

Human, Maus

Positivkontrollen:

WB : Maushirngewebe,

IP : K-562-Zellen,

Hintergrundinformationen

ZMYM3, also named as DXS6673E, KIAA0385, and ZNF261, is a 1370 amino acid protein, which may be a component of a BHC histone deacetylase complex. ZMYM3 is most abundant in brain, moderate in muscle and heart, low in other tissues except placenta. ZMYM3 plays a role in the regulation of cell morphology and cytoskeletal organization.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



mouse brain tissue were subjected to SDS PAGE followed by western blot with 25742-1-AP (ZMYM3 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.



IP Result of anti-ZMYM3 (IP:25742-1-AP, 4ug; Detection:25742-1-AP 1:300) with K-562 cells lysate 3600ug.